

DEPARTMENT OF
OPHTHALMOLOGY

Visual Field Test



This test is necessary for the eye doctor to diagnose or monitor patients with glaucoma. It usually takes about 6 to 8 minutes per eye.

The Day Before Your Test

1. Get a good night sleep!

You need rest in order to be able to concentrate and perform the test properly, so that the results are accurate and reliable.

On The Day Of Your Test

2. Instruction for test

We will provide instructions and a demo if it is your first time, or if you have forgotten how to do the test. If you have any language barriers in understanding, please ask for an interpreter.



3. Bring your current distance glasses or contact lenses

Please avoid wearing contact lenses on the day of the test. If you have to, please bring along a casing and solution for easy removal of your contact lenses during eye examination.



4. If you had a cataract operation done recently...

Inform the technician so that he or she can choose the appropriate corrective lens for your test.



5. If you have dry eyes...

We recommend that you bring your regular eye drops and apply them 10 minutes prior to test.

During The Test



6. Correct sitting position

Rest your lower back on the back rest of the chair and lean your head slightly forward. You should look straight ahead without tilting or turning your head.



7. Correct eye position

Your eye lashes should be nearly touching the focusing lens in front of your eyes. If the lens holder is blocking your view, please inform the technician.



8. Focusing on the target

Keep looking straight ahead at the target without moving your eyes. You should use your side vision to spot the lights that come into your view.


Important: The test will NOT be accurate if you move your eyes!



9. Clicking the button

Click the button only if you see a spot of light. These spots have different degrees of brightness and the time gap between each spot also changes.

Example



Click button when you see the spot of light using your side vision.

10. When you are feeling tired

You can blink or rest during the test. Press and hold the button to pause the test. The test will resume when you release the button.

1. Press and hold button down and rest your eyes quietly



Press and hold button

2. Release button and continue the test when you are ready



Release button

Clinic Appointments and Eye Screening

Tel: (65) 6357 7000

Email: contact@ttsh.com.sg

Website: www.ttsh.com.sg

LASIK Enquiries

Tel: (65) 6357 2255

Email: lasik@ttsh.com.sg

Website: www.ttshlasik.com.sg



Scan the QR Code with your smart phone to access the information online or visit <http://bit.ly/TTSHConditions-Treatments>.



Disclaimer

This is a general guide. If in doubt, please consult your doctor. Information is subject to revision without notice. The contents of this leaflet are not to be reproduced in any form without the prior permission of NHG Eye Institute.

Printed April 2023

PECC-EYE-ED-2023-1672-v1